

Microscope Compatibility for Merlin and ROLS		
<u>Brand</u>	<u>Model</u>	<u>ROLS[∞] / VROLS</u>
Alcon	LuxOR (Illumin-i)	manual (11313), auto (11363), VROLS-Z
Inami	L-0980A	manual (11313), auto (11363), VROLS-Z
Leica	M500	manual (11315), auto (11364), VROLS-L
Leica	M501	manual (11315), auto (11364), VROLS-L
Leica	M520	manual (11315), manual X (11374), auto (11364), VROLS-L, VROLS-LX
Leica	M620	manual (11315), auto (11364), VROLS-L
Leica	M650	manual (11315), auto (11364), VROLS-L
Leica	M651	manual (11315), auto (11364), VROLS-L
Leica	M655	manual (11315), auto (11364), VROLS-L
Leica	M690	manual (11315), auto (11364), VROLS-L
Leica	M691	manual (11315), auto (11364), VROLS-L
Leica	M695	manual (11315), auto (11364), VROLS-L
Leica	M820	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M822	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M840	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M841	manual (11315), auto (11364), VROLS-L; mount above beam splitter

Leica	M844	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Möller-Wedel (Haag-Streit)	Hi-R 900	11313 (manual), 11363 (auto), VROLS-Z
Möller-Wedel (Haag-Streit)	Hi-R 900 SL	11313 (manual), 11363 (auto), VROLS-Z
Möller-Wedel (Haag-Streit)	Hi-R 900 SLA	11313 (manual), 11363 (auto), VROLS-Z
Möller-Wedel (Haag-Streit)	EOS 900	11313 (manual), 11363 (auto), VROLS-Z
Möller-Wedel (Haag-Streit)	Ophthalmic 900	11313 (manual), 11363 (auto), VROLS-Z
Möller-Wedel (Haag-Streit)	Spectra 900	11313 (manual), 11363 (auto), VROLS-Z
Prescott	OMNI-FLEX	11313 (manual), 11363 (auto), VROLS-Z
Storz	Urban	DOES NOT MOUNT
Takagi	OM-8	11313 (manual), 11363 (auto), VROLS-Z
Takagi	OM-18	11313 (manual), 11363 (auto), VROLS-Z
Topcon	OMS 600	11313 (manual), 11363 (auto), VROLS-Z
Topcon	OMS 800	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	Lumera	11313 (manual), 11363 (auto), VROLS-Z; mount above beam splitter
Zeiss	MDO	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	MDU	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	MDI	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI 1	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI 6	11313 (manual), 11363 (auto), VROLS-Z

Zeiss	OPMI 6-CFC	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI 6S	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI 6-SFC	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI 6-SFR	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI 6-CH	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI CS shell (various modules below can snap in)	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	Retrolux CS module on OPMI CS base	see OPMI CS shell information
Zeiss	Retrolux 2 module on OPMI CS base	see OPMI CS shell information
Zeiss	Retrolux 3 module on OPMI CS base	see OPMI CS shell information
Zeiss	Retroskop on OPMI CS base	see OPMI CS shell information
Zeiss	OPMI Pico	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI VISU 140	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI VISU 150	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI VISU 160	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI VISU 200	11313 (manual), 11363 (auto), VROLS-Z
Zeiss	OPMI VISU 210	11313 (manual), 11363 (auto), VROLS-Z