



Optimetrics

www.optimetrics.com

sales@optimetrics.com

Toll Free: 1-866-348-3009

Tel: (305) 863-41-44



Envi S2-AR



The Envi Family of Ophthalmic Exam Chairs

Crafted with the eye care professional in mind, Envi chairs embody a perfect blend of **quality, affordability, and functionality**. Designed to meet the industry's highest standards, each chair incorporates essential features that enhance patient comfort and clinical efficiency.

- **Complete automatic power:** Effortlessly adjust chair up/down, recline back/forward, including a fully flat position at the touch of a button, and activate auto return mode.
- **Customizable comfort:** High-quality, durable upholstery is available in a variety of colors to complement your space.
- **Practical design:** Enjoy the convenience of a folding footrest with spring assist, adjustable headrest, and armrests that move synchronously with the chair.
- **Enhanced mobility:** Rotate chair 180° with the foot lock for easy positioning.
- **User-friendly controls:** Easily operate chair with soft-touch buttons conveniently located on both sides.



Technical Specifications:

- Chair base: 21" x 23.75"
- Height in lowest position: 23.6"
- Height in highest position: 31.5"
- Length when fully reclined: 65"
- Lifting capacity: 500 lbs



www.optimetrix.com
sales@optimetrix.com
Toll Free: 1-866-348-3009
Tel: (305) 863-41-44

