



Optimetrics

www.optimetrics.com

sales@optimetrics.com

Toll Free: 1-866-348-3009

Tel: (305) 863-41-44



Envi S2-MPC



The Envi Family of Ophthalmic Exam Chairs

Crafted with the eye care professional in mind, Envi chairs embody a perfect blend of **quality, affordability, and functionality**. Designed to meet the industry's highest standards, each chair incorporates essential features that enhance patient comfort and clinical efficiency.

The Envi S1-MPC is a versatile, electrically powered chair designed for a variety of medical settings. Ideal for screening rooms, laser procedures, and adaptable as a compact exam lane when paired with the Envi Switch, this chair offers exceptional functionality and patient comfort.

Equipped with an **electric power base and convenient foot pedal**, the S1-MPC effortlessly adjusts to desired heights. Patient well-being is prioritized through the inclusion of supportive handrails and a comfortable footrest, enhancing overall procedural experience. The chair's spacious upholstered backrest can be adjusted vertically to accommodate diverse patient needs and preferences.



Technical Specifications:

- Chair base: 21" x 23.75"
- Height in lowest position: 20"
- Height in highest position: 28"
- Lifting capacity: 500 lbs
- Voltage: 120v



www.optimetrics.com
sales@optimetrics.com
Toll Free: 1-866-348-3009
Tel: (305) 863-41-44

