



 **Optimetrics**

www.optimetrics.com

sales@optimetrics.com

Toll Free: 1-866-348-3009

Tel: (305) 863-41-44

Envi S2-S



Envi S2-S

The S2 Stand is a testament to our commitment to innovation and customer satisfaction. Building upon the success of the S1, we've incorporated industry-leading features and refined the design to deliver a superior instrument stand.

Experience the difference that comes from continuous improvement and a dedication to quality. Designed to pair flawlessly with the Envi line of chairs, the Series 2-S is the ultimate solution to serve the needs of the eye care professional. With easy access to controls and ring light charging wells on the control panel, the Series 2-S is the ultimate solution to serve the needs of the eye care professional.

Features:

- **Touch-Activated Controls:** Experience seamless control with the S2's touch-activated interface, enabling effortless adjustments for Envi Chair height and engaging the convenient auto-return feature.
- **Enhanced Ergonomics:** The S2 Stand prioritizes clinician well-being with its ergonomic design, minimizing strain and maximizing workflow.
- **Real-Time Instrument Readiness:** Ring light indicators provide real-time charging status for handheld instruments, ensuring they are always ready for use.
- **Compact Footprint:** The S2's compact footprint optimizes space within your exam lanes, allowing for greater ease of movement and increased efficiency.

Technical Specifications:

- Base size: 23.6" x 20.8"
- Total height of unit: 84.6"
- Slit lamp arm vertical range: 12"
- Slit lamp arm weight capacity: 80 lbs
- Phoropter arm vertical range: 10"
- Voltage: 120v

