



 **Optimetrix**

www.optimetrix.com

sales@optimetrix.com

Toll Free: 1-866-348-3009

Tel: (305) 863-41-44



Envi S2-TC



The Envi Family of Ophthalmic Exam Chairs

Crafted with the eye care professional in mind, Envi chairs embody a perfect blend of **quality, affordability, and functionality**. Designed to meet the industry's highest standards, each chair incorporates essential features that enhance patient comfort and clinical efficiency.

- **Complete power:** Easily adjust chair up/down with soft-touch buttons located on the side.
- **Customizable comfort:** High-quality, durable upholstery is available in a variety of colors to complement your space.
- **Practical design:** Enjoy the convenience of a folding footrest with spring assist and an adjustable headrest.
- **Enhanced mobility:** Rotate chair 180° with the foot lock for versatile positioning.
- **Effortless Tilt:** Control the chairs tilt function with a conveniently located handle on back of the chair.



Technical Specifications:

- Chair base: 21" x 23.75"
- Height in lowest position: 23.6"
- Height in highest position: 31.5"
- Length when fully reclined: 65"
- Lifting capacity: 500 lbs
- Voltage: 120v



www.optimetrics.com
sales@optimetrics.com
Toll Free: 1-866-348-3009
Tel: (305) 863-41-44

