




Far Vision Device

24" Wide LCD Screen, High Technology, Integrated System



- | | |
|---|-----------------------------------|
| ✓ High Contrast LCD | ✓ Wall mounting and Portable |
| ✓ Calibrated LCD | ✓ Red-Green Stereo Test |
| ✓ Full HD AntiGlare LCD 24" | ✓ Fan-Less Silent System |
| ✓ Software and hardware engineered by DMD | ✓ Free software updates unlimited |
| ✓ Settings fully customizable | ✓ Driving licence procedure |


Charts, ETDRS, Ishihara, Bio/Bino tests, Stereo Test (RG) and much more...



Specifications

- | | |
|--|---|
| ✓ TFT LCD 24" Full HD, 250 cd/m ² , 1000:1 Contrast | ✓ VA range: from 0.03 to 2.0 (0.3/10 to 20/10) |
| ✓ Patient-Screen distance: up to 8 m, step 5 cm | ✓ Unlimited number of charts, full optotype randomization |
| ✓ Maddox LED Test | ✓ Video Player Full HD with stereo audio |
| ✓ Free internet software updates, forever | ✓ Standard wall mounting Vesa 100 |
| ✓ Classification: CE, Medical Device Classe 1 | ✓ 6 Years Warranty (2 Years for LCD) |
| ✓ Compliance with: international standard UNI EN ISO 15004, ISO 10938, ISO 8596 | ✓ Dimension (WxHxD): 59.5 x 37.0 x 5.0 [cm], Weight: 4.3 Kg |

Most Useful Accessories



EyeTablet

This 10" lcd control console will provide you a perfect control of your LCD Chart Panel; no more need to look directly to the chart panel.



Light Kit

The Light Kit for Driving License test consists of two arms (right and left) on which are housed the lamps to create dazzle during testing. This is an optimized kit for the "Driving Licence Test" software supplied by us.



BlueKit

Wireless Connection Kit to Phoropters